

### TECHNICAL DATA



**R7950/5L**

#### Model R7950 Ultrasonic Couplant Gel, 250mL

- 8.45 fl oz (250ml) applicator top bottle
- Slow drying with good transducer lubrication
- Water soluble
- Includes applicator bottle
- Material Safety Datasheet (MSDS) available upon request

#### Model R7950/5L Ultrasonic Couplant Gel, 5L

- 5 Litre (1.32 gallon) barrel
- Slow drying with good transducer lubrication
- Water soluble
- Includes applicator bottle
- Material Safety Datasheet (MSDS) available upon request

#### Model R7950/12 Ultrasonic Couplant Gel, Pack of 12

- Slow drying with good transducer lubrication
- Water soluble
- Includes applicator bottle
- Material Safety Datasheet (MSDS) available upon request



**R7950**



**R7950/12**

Model	Description
R7950	Ultrasonic Couplant Gel, 250mL
R7950/5L	Ultrasonic Couplant Gel, 5L
R7950/12	Ultrasonic Couplant Gel, Pack of 12